


|  |  |
|--|--|
|  <b>Bi-Lab</b><br>QUALITY MANAGEMENT SYSTEM<br>www.rudaslab.com | ΧΗΜΕΙΟ ΠΟΥΔΑ<br>EN ISO 15189:2012 § 5.5    |
|  | Redacted by: (T.M.)<br>Approved by: (Q.M.) |

## Vitamins Profile

| Test Name  | Schedule   | Sample*             | Amount |
|--|------------|---------------------|--------|
| <b>Vitamin A</b>   | Five days  | Serum frozen        | 1.0 mL |
| <b>Vitamin B<sub>1</sub></b>                                 | Three days | EDTA Blood 4°C      | 1.0 mL |
| <b>Vitamin B<sub>2</sub></b>                                 | One week   | EDTA Blood 4°C      | 1.0 mL |
| <b>Vitamin B<sub>3</sub></b>                                 | One week   | Serum frozen        | 1.0 mL |
| <b>Vitamin B<sub>6</sub></b>                                 | One week   | Heparin/EDTA frozen | 1.0 mL |
| <b>Vitamin B<sub>9</sub></b>                                 | Daily      | Serum 4°C           | 1.0 mL |
| <b>Vitamin B<sub>12</sub></b>                                | Daily      | Serum 4°C           | 1.0 mL |
| <b>Vitamin C</b>   | One week   | Heparin frozen      | 1.0 mL |
| <b>Vitamin D (25-OH)</b>                                     | Daily      | Serum frozen        | 1.0 mL |
| <b>Vitamin D (1,25OH)</b>                                    | Five days  | Serum frozen        | 1.0 mL |
| <b>Vitamin E</b>   | One week   | Serum frozen        | 1.0 mL |
| <b>Vit B<sub>1</sub>+Vit B<sub>2</sub>+Vit B<sub>6</sub></b> | One week   | EDTA Blood + Plasma | 3.0 mL |

\* All vitamins samples must be light protected